

Resolon® Polyamid

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no. metric USP		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
1/2-Kreis schneidende = HS, mit PREMIUM-Schliff = HSM 1/2 circle, reverse cutting = HS, PREMIUM cut= HSM						
▼ HS 12 	1	5-0	0,45	blau/blue	88111	3
▽ HSM 17  PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,45	blau/blue	88121	3
▼ HS 18 	1 1.5	5-0 4-0	0,45 0,45	blau/blue blau/blue	88112 88113	3 3
▼ HS 22 	2	3-0	0,45	blau/blue	88114	3
3/8-Kreis schneidende = DS, mit PREMIUM-Schliff = DSM, DSM F, DSMF, DSMF X 3/8 circle, reverse cutting = DS, PREMIUM cut = DSM, DSM F, DSMF, DSMFX						
▽ DSM 7  PREMIUM-Schliff / PREMIUM cut	0.5	7-0	0,45	blau/blue	881543	3
▽ DSM 11  PREMIUM-Schliff / PREMIUM cut	0.7 1	6-0 5-0	0,45 0,45	blau/blue blau/blue	88151 88152	3 3
▽ DSM 11 F  PREMIUM-Schliff / PREMIUM cut	0.5	7-0	0,45	blau/blue	881528	3
▽ DSMF 11  PREMIUM-Schliff „fein“ / PREMIUM cut „fine“	0.7 1	6-0 5-0	0,45 0,45	blau/blue blau/blue	88221 882216	3 3
▼ DS 13 	0.7 1 1.5	6-0 5-0 4-0	0,45 0,45 0,45	blau/blue blau/blue blau/blue	881430 881431 881432	3 3 3

Resolon® Polyamid

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
▽ DSM 13  PREMIUM-Schliff / PREMIUM cut	0.7	6-0	0,45	blau/blue	88153	3
	1	5-0	0,45	blau/blue	88154	3
	1.5	4-0	0,45	blau/blue	881523	3
▽ DSMF 13  PREMIUM-Schliff „fein“ / PREMIUM cut „fine“	0.7	6-0	0,45	blau/blue	88222	3
	1	5-0	0,45	blau/blue	88223	3
	1.5	4-0	0,45	blau/blue	882217	3
▽ DSMFX 13  PREMIUM-Schliff Spitze „fein“ / PREMIUM cut point „fine“	1	5-0	0,45	blau/blue	882221	3
	1.5	4-0	0,45	blau/blue	882228	3
▼ DS 16 	0.7	6-0	0,45	blau/blue	881440	3
	1	5-0	0,45	blau/blue	88141	3
	1.5	4-0	0,45	blau/blue	88142	3
	2	3-0	0,45	blau/blue	881433	3
	1	5-0	0,75	blau/blue	881451	3
	1.5	4-0	0,75	blau/blue	881452	3
▽ DSM 16  PREMIUM-Schliff / PREMIUM cut	0.7	6-0	0,45	blau/blue	881529	3
	1	5-0	0,45	blau/blue	88155	3
	1.5	4-0	0,45	blau/blue	88156	3
	2	3-0	0,45	blau/blue	881548	3
	1	5-0	0,75	blau/blue	881530	3
▽ DSMF 16  PREMIUM-Schliff „fein“ / PREMIUM cut „fine“	0.7	6-0	0,45	blau/blue	882210	3
	1	5-0	0,45	blau/blue	882211	3
	1.5	4-0	0,45	blau/blue	882212	3
▼ DS 18 	1	5-0	0,45	blau/blue	88143	3
	1.5	4-0	0,45	blau/blue	88144	3
	2	3-0	0,45	blau/blue	88146	3
	3	2-0	0,45	blau/blue	881454	3
	1	5-0	0,75	blau/blue	881450	3
	1.5	4-0	0,75	blau/blue	88145	3
	2	3-0	0,75	blau/blue	88147	3
	3	2-0	0,75	blau/blue	88148	3
	1.5	4-0	4 x 0,45	blau/blue	M 881421 Multipack	2
▼ DS 18 F 	1.5	4-0	0,45	blau/blue	881441	3
	1.5	4-0	0,75	blau/blue	881442	3

Resolon® Polyamid
Sterile Nadel-Faden-Kombination
Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
▽ DSM 18  PREMIUM-Schliff / PREMIUM cut	1	5-0	0,45	blau/blue	88157	3
	1.5	4-0	0,45	blau/blue	88158	3
	2	3-0	0,45	blau/blue	881510	3
	3	2-0	0,45	blau/blue	881531	3
	1.5	4-0	0,75	blau/blue	88159	3
	2	3-0	0,75	blau/blue	881511	3
	3	2-0	0,75	blau/blue	881512	3
▽ DSM 18 F  PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,45	blau/blue	881541	3
	1.5	4-0	0,75	blau/blue	881542	3
	1.5	4-0	4 x 0,45	blau/blue	M 881548 Multipack	2
▽ DSMF 18  PREMIUM-Schliff „fein“ / PREMIUM cut „fine“	0.7	6-0	0,45	blau/blue	88224	3
	1	5-0	0,45	blau/blue	88225	3
	1.5	4-0	0,45	blau/blue	88226	3
	2	3-0	0,45	blau/blue	88229	3
	1.5	4-0	0,75	blau/blue	882218	3
	2	3-0	0,75	blau/blue	882219	3
▽ DSMFX 18  PREMIUM-Schliff Spitze „fein“ / PREMIUM cut point „fine“	1	5-0	0,45	blau/blue	882227	3
	1.5	4-0	0,45	blau/blue	882222	3
▼ DS 21 	1.5	4-0	0,45	blau/blue	88149	3
	2	3-0	0,45	blau/blue	881411	3
	1.5	4-0	0,75	blau/blue	881410	3
	2	3-0	0,75	blau/blue	881412	3
	3	2-0	0,75	blau/blue	881413	3
▽ DSM 21  PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,45	blau/blue	881513	3
	2	3-0	0,45	blau/blue	881515	3
	1.5	4-0	0,75	blau/blue	881514	3
	2	3-0	0,75	blau/blue	881516	3
	3	2-0	0,75	blau/blue	881517	3
▼ DS 24 	1	5-0	0,45	blau/blue	881458	3
	1.5	4-0	0,45	blau/blue	881424	3
	2	3-0	0,45	blau/blue	881414	3
	3	2-0	0,45	blau/blue	881425	3
	1.5	4-0	0,75	blau/blue	881457	3
	2	3-0	0,75	blau/blue	881415	3
	3	2-0	0,75	blau/blue	881416	3

Resolon® Polyamid

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
▽ DSM 24 	1.5	4-0	0,45	blau/blue	881544	3
	2	3-0	0,45	blau/blue	881518	3
	3	2-0	0,45	blau/blue	881524	3
PREMIUM-Schliff / <i>PREMIUM cut</i>	1.5	4-0	0,75	blau/blue	881545	3
	2	3-0	0,75	blau/blue	881519	3
	3	2-0	0,75	blau/blue	881520	3
▽ DSMF 24 	1.5	4-0	0,45	blau/blue	882213	3
	2	3-0	0,45	blau/blue	882214	3
	3	2-0	0,45	blau/blue	882215	3
PREMIUM-Schliff „fein“ / <i>PREMIUM cut „fine“</i>	2	3-0	0,75	blau/blue	882220	3
▽ DSMFX 24 	2	3-0	0,45	blau/blue	882223	3
PREMIUM-Schliff Spitze „fein“ / <i>PREMIUM cut point „fine“</i>						
▼ DS 27 	2	3-0	0,45	blau/blue	881453	3
	3	2-0	0,45	blau/blue	881423	3
	3.5	0	0,45	blau/blue	881456	3
	2	3-0	0,75	blau/blue	881417	3
	3	2-0	0,75	blau/blue	881418	3
▽ DSM 27 	3	2-0	0,45	blau/blue	881525	3
	3	2-0	0,75	blau/blue	881526	3
PREMIUM-Schliff / <i>PREMIUM cut</i>						
▼ DS 30 	1.5	4-0	0,45	blau/blue	881434	3
	2	3-0	0,45	blau/blue	881426	3
	3	2-0	0,45	blau/blue	881428	3
	3.5	0	0,45	blau/blue	881429	3
	1.5	4-0	0,75	blau/blue	881435	3
	2	3-0	0,75	blau/blue	881427	3
	3	2-0	0,75	blau/blue	881419	3
	3.5	0	0,75	blau/blue	881420	3
▽ DSM 30 	2	3-0	0,45	blau/blue	881532	3
	3	2-0	0,45	blau/blue	881534	3
	3.5	0	0,45	blau/blue	881536	3
PREMIUM-Schliff / <i>PREMIUM cut</i>	2	3-0	0,75	blau/blue	881533	3
	3	2-0	0,75	blau/blue	881535	3
	3.5	0	0,75	blau/blue	881537	3
▽ DSMFX 30 	2	3-0	0,75	blau/blue	882224	3
	3	2-0	0,75	blau/blue	882225	3
PREMIUM-Schliff Spitze „fein“ / <i>PREMIUM cut point „fine“</i>						

Resolon® Polyamid
Sterile Nadel-Faden-Kombination
Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
▼ DS 40 	3	2-0	0,45	blau/blue	881436	3
	3.5	0	0,45	blau/blue	881438	3
	3	3-0	0,75	blau/blue	881437	3
	3.5	0	0,75	blau/blue	881439	3
	3	2-0	0,90	blau/blue	881455	3
▽ DSM 40 	3	2-0	0,75	blau/blue	881538	3
	3.5	0	0,75	blau/blue	881540	3
PREMIUM-Schliff / <i>PREMIUM cut</i>						
▽ DSMFX 40 	3	2-0	0,75	blau/blue	882226	3
PREMIUM-Schliff Spitze „fein“ / <i>PREMIUM cut point „fine“</i>						
Gerade Rundkörper, biegsame Nadel = GRW <i>Straight round bodied, flexible needle = GRW</i>						
● 2 x GRW 65 	3	2-0	0,75	blau/blue	88091	3
	3.5	0	0,75	blau/blue	88092	3
Biegsame Nadel / <i>Flexible needle</i>						
Gerade schneidende = GS <i>straight, cutting = GS</i>						
▼ GS 60 	2	3-0	0,75	blau/blue	88171	3
	3	2-0	0,75	blau/blue	88172	3

Resolon® Polyamid

Einzelfadenpackungen, steril in vorgeschnittenen Längen

Sterile pre-cut suture, non-needed

Stärke Nr. gauge no. metric USP		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
1.5	4-0	3 x 0,50	blau/blue	X 335 N	2
2	3-0	3 x 0,50	blau/blue	X 435 N	2
3	2-0	3 x 0,50	blau/blue	X 535 N	2
3.5	0	3 x 0,50	blau/blue	X 635 N	2
2	3-0	10 x 0,50	blau/blue	X 4105 N	2